PRIMAX[®] SB-X

Fluoroscopic medical X-ray film for mass chest X-ray diagnostics

Applications

Orthochromatically sensitized film, suitable for use in fluoroscopic cameras (for instance Carl Zeiss RK-70 or RK-11, Russian cameras KF 70 or KF 400), working with green-emitting screens (for instance CAWO LS 16, RENEX ERS-S-2A or Renex EFG-G-2B).

Density curve

see diagram

Processing

a) Machine processing at 33-35°C

b) Manual processing at 21-24°C

It is to be considered that intermediate rinsing will be required for 30 seconds after developing and for 30 minutes after fixing to ensure longest possible storage of the processed films.

Film properties

- excellent detail sharpness
- optimal contrast range
- low patient dose due to high speed

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
Remark

The ratio between temperature of the developer and developing time as shown in the right diagram above has to be considered. Any ratio within the marked field may be selected. The higher the temperature of the developer the shorter the developing time.

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Mammography Film PRIMAX X-MA

Applications:

The high speed and high contrast medical X-ray film type PRIMAX X-MA is specially designed for routine mammography examinations

Characteristics/Handling/Processing:

- Single-coated orthochromatic film on blue tinted PET-base
- For easy handling of the film in the darkroom the film sheets are fitted with a notch. Viewing at the coated film side the notch is at the upper right edge. The coated film side must be directed towards the intensifying screen during the exposure of the film
- Automatic processing in standard 90 second cycle in any type of standard branded X-ray film processing chemicals

Intensifying screens/Packing:

- green emitting mammography screens to be used
- packing in boxes of 100 sheets NIF in sizes 18x24 and 24x30 cm

Storage:

- unprocessed film should be stored at maximum 23°C
- relative humidity to be kept within 40 to 60 per cent
- protection necessary against ionizing radiation (max. 90 nGy/h)
- protection necessary against the influence of chemicals

PRIMAX[®] K-RH/K-RI

Laser Films PRIMAX K-RH and K-RI

Applications:

For high quality diagnostic recordings from laser printers. Each type is matching with its optimum sensitivity to the emitted exposure light of the corresponding laser printer

Characteristics/Handling/Processing:

- Single emulsion films
- type K-RH for printers with sensitivity to red light at 633 nm, type K-RI with sensitivity to infrared light at 820 nm
- Generally for darkroom loading, but also suitable for daylight loading in selected Kodak printers
- Automatic processing in standard 90 second cycle in any type of standard branded X-ray film processing chemicals. Processing possible also in combination with our standard double coated screen films

Packing and shelf life:

- packing in boxes of 100 sheets NIF in size 14x17" (35x43 cm)
- shelf life 18 months for type PRIMAX K-RI, 36 months for type PRIMAX K-RH

Storage:

- unprocessed film should be stored at maximum 23°C
- relative humidity to be kept within 40 to 60 per cent
- protection necessary against ionizing radiation (max. 90 nGy/h)
- protection necessary against the influence of chemicals

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