



MicroDrop® Family2

Aerosol therapy for the whole family



MicroDrop® Family2 is a compressor operated nebulizer system for the treatment of respiratory diseases of the lower airways. Combined with the Baby-Kit the nebulizer can already be used for infants.

The **RF6 plus nebulizer** can be used for all approved and tested inhalation solutions and suspensions.

- *Short treatment time: 2.5 ml in 8 minutes*
- *Valve system for mouthpiece- and mask-inhalation*
- *Anti-spilling system – no medication can spill out*
- *Compressor for continuous use*

Delivery includes: Compressor, RF6 plus nebulizer, mouthpiece, children and adult mask, tubing, compressor filter, carrying bag, manual



Aerosol characteristics*

Particle size	3.34 µm (MMAD)
Total output rate	320 mg/min
Breathable fraction < 5 µm	69.2 %
RDDR < 5 µm	44.22 µg/min (open valve) 46.36 µg/min (closed valve)
RDDR < 3 µm	27.18 µg/min (open valve) 30.64 µg/min (closed valve)



Technical data

Pressure compressor	approx. 0.6 bar
Weight	1.2 kg
Size (H x L x W)	105 x 200 x 170 mm
Noise level	54 dB(A)
Filling volume nebulizer	min. 2 ml / max. 8 ml
Compressor air output	appx. 9 L/min

DESCRIPTION

REF

MicroDrop® Family2

M51604-00

* In vitro testing by TÜV Rheinland LGA Products GmbH – Germany in compliance with European Standard EN 13544-1 for nebulizing systems.