

**X-RAY PROTECTION**

**PATIENT'S PROTECTION**

**PRIMAX<sup>®</sup> radiation-protective clothing  
for use in medical radiology**



**Gonad protection  
PRIMAX<sup>®</sup> soft**

(Half apron, 0,50 mm Pb)

- for protection of patient's gonads during radioscopy of thorax and upper abdomen
- available in 4 sizes: 25x20, 30x30, 37x45 and 40x50 cm

**Heavy protective neck cover PRIMAX<sup>®</sup> soft  
with extension**

(0,50 mm Pb, in addition to protection aprons)

- suitable for being worn in combination with protective PRIMAX Soft aprons
- used for protection of thyroid and sternum
- easy to handle neck Velcro-fastener
- neck cover suitable for size 38-46, collar height 6 cm



**Packing, transport, storage, maintenance  
and disposal**

- detailed information is given in our directions for use which are attached to each article
- all PRIMAX articles are individually packed in stable paper board containers
- storing temperature should be kept between - 5 and + 35°C
- all articles are easily washable and sanitisable
- regular inspections for damages before and regular cleaning after use is necessary
- waste disposal of unusable articles to made on the base of key number 170403 European Waste Catalogue

**Quality assurance**

- PRIMAX Berlin GmbH is maintaining a quality management system as per ISO 13485:2003
- PRIMAX soft radiation-protective articles comply with all guidelines of IEC 61331-3:1998 as well as directive 89/686/EEC (articles for protection of medical personnel) and directive 93/42/EEC (articles for patient protection) and are therefore labelled with the CE mark



**CE**  
0123



### **DENTAL APRON WITH INTEGRATED THYROID PROTECTOR**

(0,35 mm Pb)

- full coverage of front body and shoulder of the patient
- Velcro-fastener in the neck (see above)
- size: 60 x 80 cm



### **DENTAL APRON WITHOUT THYROID PROTECTOR**

(0,35 mm Pb)

- full coverage of front body and shoulder of the patient
- Velcro-fastener in the neck (see above)
- size: 60 x 80 cm



### **DENTAL APRON WITH VERTEBRAL PROTECTION**

(0,35 mm Pb)

- for panorama and longdistance X-rays
- special protection of the backbone
- convenient Velcro-fastener on the shoulder
- size: 60 x 80 cm